

Long Exposure

Gene Koshinski and Tim Broschious

♩ ≈ 110

Switch held in left hand, pointing upward
Right hand begins resting on leg
Head stationary (looking forward)

A

Gene

Switch held in left hand, pointing upward
Right hand begins resting on leg
Head stationary (looking forward)

Tim

Sample

B

Gene

*Follow the tip of the switch with your head quickly down, then with your head/eyes as you slowly move the switch back to resting position (then back to looking forward)

Tim

*Follow the tip of the switch with your head quickly down, then with your head/eyes as you slowly move the switch back to resting position (then back to looking forward)